

I-CONNECT STUDENT

Procedural Fidelity

Student ID: _____

Mentor: _____

Observer Name: _____

Observer2/Reliability: _____

Date: _____

Start Time End Time TotalTime____

I-CONNECT Procedures	Observed	Quality
1) Student unlocked device	Y N	1 2 3
2) Student opened I-Connect	Y N	1 2 3
3) Student started I-Connect		
a) Logged in	Y N	1 2 3
b) Selected 'My Monitor'	Y N	1 2 3
c) Selected 'Location'	Y N	1 2 3
d) Selected 'Start Monitoring'	Y N	1 2 3
e) Selected 'Pause' & 'Resume' Monitoring	N/A Y N	1 2 3
f) Selected 'Stop Monitoring' at end of session	Y N	1 2 3

Total Fidelity Score_____

Total Quality Score_____

Total Score Possible_____

Total Score Possible_____

Total Score divided by Total Possible = %yes_____

Average_____

If fidelity does not reach 80% booster session indicated Booster session indicated Yes No

If fidelity does not reach 80% booster session indicated

Booster session indicated Yes No